

Abaxis Driver CD Version 1.3.0 Release Notes

About

This software release note is applicable to the Abaxis Driver CD Version 1.3.0.

The Abaxis Driver CD contains Windows, Macintosh and Linux Virtual Communications Port drivers for the Piccolo Xpress and VetScan VS2 analyzers.

Virtual Communications Port drivers found on the Abaxis Driver CD are designed to enable communication between Piccolo Xpress and VetScan VS2 analyzers and Windows, Macintosh and Linux computers. Piccolo Xpress and VetScan VS2 analyzers are installed as a Virtual Communications Port on the computer.

New Features

The Abaxis Driver CD contains a 64-bit Virtual Communications Port driver for Windows XP, Vista, 7, Server 2003, Server 2008 and Server 2008 R2.
The Abaxis Driver CD contains a Virtual Communications Port driver for Macintosh OS X 10.5 (Leopard) and 10.6 (Snow Leopard).

Known Issues

There are no known issues.